

LCD Digital Orbital Shaker SK-0330-Pro

Application

Suitable for mixing in culture dishes, flasks and beakers. Ideal for aeration and extraction applications. It can be used in environmental chambers and Co2 incubators.



Specification	SK-0330-Pro
Orbital diameter	10mm
Max. load capacity (with platform)	7.5Kg
Motor type	Brushless DC motor
Motor input power	28W
Motor output power	15W
Speed range	100 - 500rpm
Speed display	LCD
Timer display	LCD
Time setting range	1min - 99h 59min
Voltage, Frequency	100 - 240V, 50/60Hz
Power	30W
Weight	12Kg
Permissible ambient temperature and humidity	5 - 40°C, 80%RH
Dimension [W x D x H]	370 x 420 x 100mm
Protection class	Ip21
Data connector	Rs232

Features

- Independent LCD display for speed and time
- Compact footprint
- Provides smooth and uniform shaking
- Two options: 7.5Kg with platform speed range 100 - 500rpm and 100 - 800rpm
- Continuous or timer operation mode
- Maintenance - free brushless DC motor
- Over speed detection and protection
- Wide range of optional platform for flasks, and other popular vessels
- Remote function provides PC control and data transmission
- Free software for PC control & data logging



Manufactured by:







DLAB SCIENTIFIC CO.,LTD

31 Yu An Road, Beijing Airport Economic Core Zone,
Shunyi District, Beijing 101318, China

Distributed by:

PT INNOVASI DIAGNOSTIKA

Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,
Cikarang Utara 17530, Indonesia

Cat. No.	Description
 18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension [W x D]: 30 x 28.5cm
 18202835	Silicone mattress for 16 x 5 - 15mL tubes, use with SK330.5
 18202846	Silicone mattress for 6 x 50mL tubes, use with SK330.5
 18900027	SK330.1, Universal attachment with 4 bars, Size: 32 x 32cm, Capacity for flask: 25mL x 21pcs, 50mL x 15pcs, 100mL x 12pcs, 250mL x 6pcs, 500mL x 4pcs
 18900040	SK330.3, Lengthways roller attachment, Size: 32 x 32cm, Capacity for flask: 25mL x 21pcs, 50mL x 15pcs, 100mL x 12pcs, 250mL x 6pcs, 500mL x 4pcs
 18900036	Bar for SK330.1 universal attachment
 18900280	Connect sticks for double dishes shaker

Manufactured by:







DLAB SCIENTIFIC CO.,LTD

31 Yu An Road, Beijing Airport Economic Core Zone,
Shunyi District, Beijing 101318, China

Distributed by:

PT INNOVASI DIAGNOSTIKA

Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,
Cikarang Utara 17530, Indonesia

Cat. No.	Description
 18900028	SK330.2, Fixing clip attachment, Size 32 x 32cm
 18900029	Fixing clip for flask volume 25mL, used with SK330.2 [9pcs], SK180.2 [28pcs]
 18900030	Fixing clip for flask volume 50mL, used with SK330.2 [9pcs], SK180.2 [16pcs]
 18900031	Fixing clip for flask volume 100mL, used with SK330.2 [9pcs], SK180.2 [16pcs]
 18900032	Fixing clip for flask volume 200 / 250mL, used with SK330.2 [4pcs], SK180.2 [9pcs]
 18900033	Fixing clip for flask volume 500mL, used with SK330.2 [2pcs], SK180.2 [9pcs]

Manufactured by:

DLAB SCIENTIFIC CO.,LTD

31 Yu An Road, Beijing Airport Economic Core Zone,
Shunyi District, Beijing 101318, China

Distributed by:

PT INNOVASI DIAGNOSTIKA

Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,
Cikarang Utara 17530, Indonesia