

Palm Micro Centrifuge D1008

Application

This centrifuges is ideal for quick spins, microfiltration of samples, cell separations and many other routine laboratory procedures.

Features

- Easy to operate with dual start/stop function close lid for a spin speed up to 7000 rpm
- Smooth and quiet operation (noise level $\leq 45\text{dB}$)
- Accelerates and brakes in seconds
- Safe, durable and robust rotor with clamp locking design
- Compatible with 2mL x 8 and PCR 8 x 4 rotors
- Quick-change rotor without any tool



Quiet Operation



Vibration Free



Specification	D1008
Max. Speed	7000rpm
Max. RCF	2680 * g
Rotor Capacity	0.2/0.5/1.5/2mL * 8 0.2mL * 32PCR tubes 0.2mL * 4PCR8-strips
Run Time	Continuous Operation
Driving Motor	DC motor
Power Requirements	AC100-240V/50Hz/60Hz20W
Noise Level	$\leq 45\text{dB}$
Dimension [D x W x H]	160 x 170 x 122mm
Weight	0.5kg
Certification	CE cTüVus FCC MCA

Manufactured by:

DLAB SCIENTIFIC CO.,LTD

31 Yu An Road, Beijing Airport Economic Core Zone,
Shunyi District, Beijing 101318, China

Distributed by:

PT INNOVASI DIAGNOSTIKA

Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,
Cikarang Utara 17530, Indonesia

Available Rotor



A4-PCR8P

Cat. No. 19201165

Max. speed: 7000rpm
 Max. RCF: 2680 x g
 Rotor capacity: 0.2mL x 32PCR tubes,
 0.2mL x 4PCR8-strips
 Rotor material: High strength plastic



A8-2P

Cat. No. 19201164

Max. speed: 7000rpm
 Max. RCF: 2680 x g
 Rotor capacity: 2mL/1.5mL x 8, 0.2mL
 x 8, 0.5mL x 8
 Rotor material: High strength plastic
 Available rotor adapters: SA02P2, SA05F



Available Adapters

SA02P2

Cat. No. 19400029

0.2mL rotor adapter
 Used with A8 - 2P rotor
 8pcs/pk



SA05P2

Cat. No. 19400026

0.5mL rotor adapter
 Used with A8-2P rotor
 8pcs/pk



Manufactured by:
DLAB SCIENTIFIC CO.,LTD
 31 Yu An Road, Beijing Airport Economic Core Zone,
 Shunyi District, Beijing 101318, China

Distributed by:
PT INNOVASI DIAGNOSTIKA
 Kawasan Industri Jababeka III, Jl. Tekno 1 Blok C2 Unit C,
 Cikarang Utara 17530, Indonesia